

© 2010 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE representative for the most current information.

GE, GE Monogram, and DPX NT are registered trademarks of General Electric Company.

Windows is a registered trademark of Microsoft Corporation.

GE Medical Systems, Ultrasound & Primary Care Diagnostics, LLC, a General Electric company, doing business as GE Healthcare.

Indications for use: The DPX Series NT Bone Densitometer provides an estimate of BMD at the spine, proximal femur and forearm regions. This BMD value can then be compared to a reference population at the sole discretion of the physician.

CAUTION: Federal Law restricts this device to sale by or on the order of a physician.

Lunar Product Division Americas

GE Healthcare Lunar
Global Headquarters
PO Box 7550
Madison, WI 53707-7550
T: +1-888-795-8627 (option 9, then option 9)
F: +1-608-223-2482

Lunar Product Division Europe, Middle-East & Africa

GE Healthcare Lunar
Kouterveldstraat 20
B-1831 Diegem, Belgium
T: +32-2-7197217
F: +32-2-7197359
info.lunar@ge.com

Lunar Densitometry Asia & Pacific Headquarters

GE Healthcare Lunar
3/F GE China Technology Park
No. 1 Hua Tuo Road
Shanghai 201203, China
T: +86-21-38777888 (Ext. 60128 or 60480)
F: +86-21-38777451

About GE Healthcare

GE Healthcare provides transformational medical technologies and services that are shaping a new age of patient care. Our broad expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, biopharmaceutical manufacturing technologies, performance improvement and performance solutions services help our customers to deliver better care to more people around the world at a lower cost. In addition, we partner with healthcare leaders, striving to leverage the global policy change necessary to implement a successful shift to sustainable healthcare systems.

Our **healthymagination** vision for the future invites the world to join us on our journey as we continuously develop innovations focused on reducing costs, increasing access and improving quality and efficiency around the world. Headquartered in the United Kingdom, GE Healthcare is a \$17 billion unit of General Electric Company (NYSE: GE). Worldwide, GE Healthcare employs more than 46,000 people committed to serving healthcare professionals and their patients in more than 100 countries. For more information about GE Healthcare, visit our website at www.gehealthcare.com

To receive **Lunar News** and be informed about the latest developments in Densitometry, **please register for SmartMail**. You will find the SmartMail registration link on the left side of all pages of www.gehealthcare.com



GE imagination at work



BMD-0179-01.10-EN-GL

GE Healthcare

Lunar DPX NT

Office-based bone densitometry



High performance bone densitometer

Bone densitometry has reached a landmark in clinical performance with the Lunar DPX NT® densitometer.


With the combination of clinical confidence and dedicated utility, easy-to-use enCORE software and high performance low dose scanning, Lunar DPX NT assists physicians to confidently and efficiently diagnose osteoporosis and assess fracture risk.

Also available for the Lunar DPX NT are innovative tools to increase productivity, such as the DualFemur feature for seamlessly scanning both femorae in one automatic process, and Composer for automated generation of custom reports.




GE Healthcare Madison, WI

Patient:	XXXXXXXX	Facility ID:	XXXXXXXX
Birth Date:	XX/XX/XXXX 41.3 years	Referring Physician:	XX/XX/XXXX
Height / Weight:	67.0 in. 137.0 lbs.	Measured:	XX/XX/XXXX 3:53:08 PM (6.00)
Sex / Ethnic:	Female White	Analyzed:	XX/XX/XXXX 3:14:19 PM (12.00)




Matched for Age, Weight (females 25-100 kg), Ethnic USA (Combined NHANES (ages 20-30) / Lunar (ages 20-40)) AP Spine Reference Population (v106)
Statistically 68% of repeat scans fall within 1SD (± 0.010 g/cm² for AP Spine L1-L4)

Image not for diagnosis



Matched for Age, Weight (females 25-100 kg), Ethnic USA (Combined NHANES (ages 20-30) / Lunar (ages 20-40)) Femur Reference Population (v110)
Statistically 68% of repeat scans fall within 1SD (± 0.012 g/cm² for Left Femur Total)



Image not for diagnosis



Matched for Age, Weight (females 25-100 kg), Ethnic USA (Combined NHANES (ages 20-30) / Lunar (ages 20-40)) Femur Reference Population (v110)
Statistically 68% of repeat scans fall within 1SD (± 0.012 g/cm² for Right Femur Total)

Image not for diagnosis

Connectivity and productivity solutions



enCORE software

enCORE software – seamless osteoporosis management


The intuitive graphical enCORE software provides a dedicated palette of clinical applications for seamless osteoporosis assessment in a streamlined operator-friendly package, while ensuring clinical confidence and fast throughput. The embedded AutoAnalysis feature facilitates the assessments even more, by automatically placing the (adjustable) standard regions of interest (ROI).



DualFemur with values per ROI

DualFemur – identifying the weakest femur


With the DualFemur option, both femorae are automatically scanned in one seamless acquisition without repositioning the patient. As such, DualFemur allows you to assess the density of the critical hip region, including identification of the weakest side increasing confidence in your treatment decisions. In addition, the trending function enables seamless follow-up of change over time.^{1,2,3}



ScanCheck

ScanCheck – add quality and diagnostic power


ScanCheck automatically studies acquisition inputs and the acquired image, looking for errors and patient irregularities. When it detects anomalies, it displays explanations and instructions which can be reviewed by the interpreting physician. ScanCheck helps speed throughput and reduces errors. It also helps technologists provide information to the interpreting physician.^{4,5}



OneScan – three sites in one test

OneScan – three sites in one test

OneScan simplifies BMD testing by acquiring lumbar spine and bilateral femur scans in one, automatic process from a single patient position, without compromising diagnostic confidence.^{6,7,8}




With SmartScan

Without SmartScan

SmartScan – increase throughput

Patient scanning speed is improved with the SmartScan feature: this time-saving scanning technique automatically adapts the scan path to the skeletal structure, imaging the necessary anatomy only, resulting in further speed gains and dose reduction. Utilizing this feature, AP spine and femur scan times are typically less than two minutes in standard mode.



Composer: custom reports

With clinical diagnosis and treatment decisions based on a variety of pre-defined criteria and guidelines established by international and local societies,¹⁰ it might not always be that easy for your referring colleagues and administrations to interpret multi-page reports. Composer allows you to automatically generate concise custom patient reports, including imagery, clinician diagnosis and monitoring assessments, in full accordance with the pre-defined criteria and guidelines in your locality.

OneVision: the spine and both femurs in one comprehensive report

OneVision automatically combines scans of the spine and both hips into one comprehensive report, acquired in one process and evaluated in one analysis. Rather than receiving multiple assessment reports, the referring physician receives a single, consolidated report that combines all risk assessment analyses.


Complete connectivity⁹ with DICOM, HL7, MUDB and TeleDensitometry

DICOM and HL7 connectivity seamlessly integrate densitometry results with Picture Archival and Communication Systems (PACS) and Radiology/Hospital Information Systems (RIS/HIS) respectively.


The Worklist feature, present in both DICOM and HL7, enables you to automatically use patient information from scheduling applications, increasing throughput while also helping to reduce data entry errors.

Multi-User Database access (MUDB) improves flexibility and productivity by offering the possibility to access and/or reanalyze scans remotely and to share with clinical partners.

TeleDensitometry enables you to send paperless reports as easy e-mail attachments, viewable on any personal computer without special software.




Dedicated to BMD and beyond



Total body BMD - body composition

Total body, body composition


Body composition measurement with dual-energy X-ray absorptiometry (DXA) can look beyond weight and the traditional body mass index (BMI) to determine body fat distribution. Body composition scans with the Lunar DPX NT provide exceptionally precise and accurate data on bone and tissue composition, including bone mineral density (BMD), lean tissue mass, and fat tissue mass. They provide both total body data and regional results (trunk, arms, legs, pelvis and android/gynoid regions). The measurements are fast and non invasive.



Pediatric

Lunar DXA pediatric application

Now you can use one powerful set of tools to get valuable clinical information about growth and development in children. The Lunar DXA pediatric application measures more than BMD - it provides assessment of bone, fat and lean tissue composition. These measurements enable enhanced evaluation of growth and development that include height for age (bone length)¹¹, BMC for bone area (bone mineralization)¹¹, bone area for height (bone width)¹¹, lean body mass for height (muscle development)^{12,13}, and BMC for lean body mass (muscle-bone balance)^{12,13}




Hip Axis Length & Femur Strength Index

Advanced Hip Assessment (AHA)

The AHA application provides tools to evaluate the structural properties of the hip:

- **Hip Axis Length (HAL)** has been demonstrated in prospective studies as an effective adjunct to femur bone density in predicting fracture risk.
- **Cross-Sectional Moment of Inertia (CSMI) and Femur Strength Index (FSI)** are calculated for the assessment of the load-bearing capacity of the hip.
- **Color bone mapping** is displayed to differentiate areas of cortical and high/low density trabecular



Orthopedic - hip implant

Orthopedic - Peri-prosthetic hip implant

The orthopedic application provides highly accurate and precise bone mineral density and bone mineral content values. Bone assessment in the vulnerable region surrounding an implant is now possible. This application also enables automated bone assessment of the hip implant using standard Gruen zones (7 zones) and extended Gruen zones (19 zones) to provide exceptional evaluation for practitioners and clinical researchers specialized in the fields of orthopedics and surgery.



Forearm

Bone evaluation of peripheral sites

The optional peripheral applications, such as the radius and ulna can be evaluated to provide additional clinical information on the skeletal status of your patient, or patient population.

Lunar DPX NT technical specifications:^{9,14}

Available applications and options

- AP spine
- Femur
- DualFemur
- OneScan
- Advanced Hip Assessment (AHA)
- ScanCheck
- Total body/body composition
- Estimated Total Body %Fat
- Forearm
- Orthopedic
- Pediatric
- OneVision
- Composer
- Lateral spine BMD
- TeleDensitometry (e-mail)
- HIPAA SecureView
- Practice Management tools
- DICOM (worklist, color print and store)
- HL7 bidirectional interface
- Multi-User Database access (MUDB) (1-3 or 1-10 users)
- SQL database
- Applaud CD-based training
- Remote connectivity for direct customer support

enCORE Windows-based user interface

- Advanced intuitive graphical interface with multimedia on-line help
- Multiple languages available
- SmartScan for scan window optimization and dose reduction
- Automated scan mode selection

- AutoAnalysis for better precision
- Customized analysis for clinical flexibility
- Exam comparison process
- Multiple patient directories with database
- BMD or sBMD results, BMC and area
- Extensive reference data: >12,000 USA/Northern European subjects, as well as NHANES, and numerous regional databases.
- T-score, Z-score, % young adult and % age matched
- WHO guidelines for diagnosis of osteoporosis
- Patient trending with previous exam importation
- enCORExpress mode for brief click path

Standard features

- Washable table pad

Quality assurance

- Automated test program with complete mechanical and electronic tests
- Automated QA trending with complete storage

Scanning method

- DXA pencil-beam technology with SmartScan technology
- No scout scan required, no moving table

X-ray characteristics

- Constant potential source at 76kV
- Dose efficient K-edge filter
- Tube current: 0.05 - 1.50 mA

Detector technology

- NaI PM tube detector
- High pulse rate

Dimensions (L x W x H) and weight

- System: 2.42m x 1.03m x 1.28m - 272kg (95" x 41" x 50" - 599lbs)
- Table height: .63m (25")

Patient weight limit

- 136kg (300lbs)

External shielding

- Not required: X-ray safety requirements may vary by location. Please inquire with local regulatory authorities.
- Operating scatter: < 0.2 mR/hr (2 µSv/hr) @ 1m (39") from X-ray source
- GE Healthcare recommends consulting your local regulatory agency to comply with local ordinances.

Environmental requirements


- Ambient temperature: 18-27°C (65-81°F)
- 120 VAC 50-60 Hz 20A dedicated circuit or 230-240 VAC 50-60Hz 10A dedicated circuit ±10%
- Humidity: 20%-80%, non-condensing

Computer workstation

- Windows platform
- Computer, printer and monitor



Minimum room dimensions:



References

1. M Kamimura, H Hirabayashi, M Konishi, Q Zhou, H Kato. Osteoporosis diagnosis and treatment decisions with Dual Femur in Japanese women. Presented at the 17th International Bone Densitometry Workshop, Kyoto Japan, November 2006
2. RE Cole. Dual Femur densitometry - Effect on diagnosis and treatment decisions. Presented at the American Society for Bone and Mineral Research Annual meeting, October 2004.
3. RE Cole. The effect of DualFemur scanning on Osteoporosis diagnosis and treatment decision making.
4. EN Schwartz, DM Steinberg. Performance of Computer Assisted Densitometry (CAD) in Spine and Femur Analysis: Comparison with Visual Assessment by Experienced Densitometrists. Abstract Published Osteoporos Int 15 (Suppl 1):S39. Poster Presented at the IOF World Congress on Osteoporosis, May 2004.
5. HS Barden, P Markwardt, R Payne, B Hawkins, M Frank, KG Faulkner (2003) Automated assessment of exclusion criteria for DXA lumbar spine scans. J Clin Densitom 6:401-409.
6. C Simonelli, L Del Rio, N Binkley. Comparison of Spine BMD Measurements from DXA With and Without Leg Elevation. Abstract Published J Bone Miner Res (2004) 19 (Suppl 1):S364. Poster Presented at ASBMR Annual Meeting, October 2004.
7. M Kamimura, H Hirabayashi, M Konishi, Q Zhou, HS Barden, H Kato. Comparison of lumbar spine BMD and T-scores with conventional and OneScan leg positioning in a Japanese population. Presented at the 17th International Bone Densitometry Workshop, Kyoto Japan, November 2006
8. RH Nord, DL Ergun, KG Faulkner. Effect of patient positioning devices on bone density measurements. Abstract Published J Bone Miner Res (2002) 17 (Suppl 1): S313. Poster Presented at ASBMR Annual Meeting, September 2002.
9. Networking is under the user's responsibility
10. The World Health Organization (WHO), the International Society of Clinical Densitometry (ISCD) and the National and International Osteoporosis Foundation (NOF and IOF)
11. Molgaard C, Thomsen BL, Prentice A, Cole TJ, Michaelsen KF (1997) Arch Dis Child 76:9-15.
12. Crabtree NJ, Kibirge MS, Fordham JN, Banks LM, Muntoni F, Chinn D, Boivin CM, Shaw NJ (2004) The relationship between lean body mass and bone mineral content in paediatric health and disease. Bone 35:965-972.
13. Schoenau E, Neu CM, Beck B, Manz F, Rauch F (2002) Bone mineral content per muscle cross-sectional area as an index of the functional muscle-bone unit. J Bone Miner Res 17:1095-1101.
14. Depending on product configuration and availability. Contact GE Healthcare or our local distributor for the detailed current configuration and optional hardware.



800-356-3388
978-374-6371
Fax - 978-521-2214

49 Newark Street
Haverhill, MA 01832
sales@associatedxray.com
www.associatedxray.com