


Membrane
Console



Compact High
Frequency
Generator

4-way Float Top Table

Hands Free Collimator Light
Control and Table Lock Release

Standard Vet



When you need the best
Analog X-ray system

- Membrane Console
- 32kW/400mA Compact HF Generator
- Hands Free Collimator light control
- Variable SID for table top and cassette tray imaging
- Rotating Tube Head and collimator
- X-ray tube Rotor Brake
- 4 Way Float Top Table - 5ft Standard
- 208/230/240 VAC Single Phase



**SIMPLE
DEPENDABLE
SMART**

Standard Vet System Dimensions



Standard X-ray Tube:

- 1.0/2.0 Focal Spots
- 16° Target angle
- 140,000 HU

Other Option Available:

- Custom (38" to 78") length top
- Fractional Focal Spot X-ray tubes
- 3 phase and battery powered generators

Grid Cabinet:

- 8:1 103 line pair grid
- Stainless steel cassette tray



MEDICAL X-RAY SYSTEMS
EQUIPMENT • SUPPLIES • SERVICE
800-356-3388

Associated X-Ray Imaging Corp.
49 Newark Street
Haverhill, MA 01832

Toll Free: 800-356-3388
Phone: 978-374-6371
Fax: 978-521-2214

www.associatedxray.com
sales@associatedxray.com

**SIMPLE
DEPENDABLE
SMART**