

Superior diagnostic imaging and procedural guidance in a wireless handheld ultrasound system.



POCKETPRO H<sup>2</sup> Handheld Ultrasound System

The POCKETPRO H<sup>2</sup> 500L linear wireless handheld ultrasound system is designed for MSK, Vascular Access, and procedural guidance.

- Range of 6 to 12 MHz
- Color, PW, M-mode functions
- Customizable presets
- JPEG and DICOM file formats

# Confident performance with every scan.

## Mobile app

Supported by the SononX app, available on Apple, Android and Windows platforms.


- Intuitive user interface supports one-hand operation and ease-of-use
- On-board patient data management solutions provide connected patient care
- Seamless connectivity with PACs

## Seamless connectivity


The handheld ultrasound device is designed for secure and seamless connectivity with PACs to complement a facility's workflow. Secure enterprise environments are not compromised since there is no requirement to store images outside your organization, no annual subscription, no device license limit and you retain full ownership of your images.

## Increased ROI


Combined with industry leading support, online education and personal one-on-one applications training, the POCKETPRO H<sup>2</sup> delivers a strong return on your investment in the wireless handheld ultrasound market.




Muscle




Tendon



Bone Cortex



Nerve



Bursa



Ligament

**BOOK A LIVE ULTRASOUND DEMO WITH AN EXPERT! SCAN THE CODE OR VISIT [calendly.com/demoskmha/ultrasound-virtual-demos](https://calendly.com/demoskmha/ultrasound-virtual-demos)**

Reliable and affordable diagnostic imaging and procedural guidance for point-of-care professionals.

- Observe live ultrasound scanning • See the difference superior image quality makes
- Learn the advantage of wireless scanners • Ask about affordability and financing

Manufactured for Konica Minolta Healthcare

© 2023 Konica Minolta Healthcare Americas, Inc. SONON is a registered trademark of Healcerion USA.



**MEDICAL X-RAY SYSTEMS**  
EQUIPMENT • SUPPLIES • SERVICE  
800-356-3388

Associated X-Ray Imaging Corp.  
49 Newark Street  
Haverhill, MA 01832  
Toll Free: 800-356-3388  
Phone: 978-374-6371  
Fax: 978-521-2214

[www.associatedxray.com](http://www.associatedxray.com)  
[sales@associatedxray.com](mailto:sales@associatedxray.com)