



# Acuity CDR



# R<sub>admedi</sub>X

# R<sub>admedi</sub>X

## CHIROPRACTIC X-RAY REIMAGINED



**MEDICAL X-RAY SYSTEMS**  
EQUIPMENT • SUPPLIES • SERVICE  
800-356-3388

Associated X-Ray Imaging Corp.  
49 Newark Street  
Haverhill, MA 01832

Toll Free: 800-356-3388  
Phone: 978-374-6371  
Fax: 978-521-2214

[www.associatedxray.com](http://www.associatedxray.com)  
[sales@associatedxray.com](mailto:sales@associatedxray.com)



## Acuity CDR, Advanced X-Ray Technology

Radmedix has expanded our proven straight arm line with a new advanced x-ray system designed specifically for chiropractic care, the Acuity CDR. The Acuity CDR is an ultra small footprint, compact system which comes with a variety of features such as full software integration, 20° arm rotation, magnetic lock brakes, 40kw generator, LED collimator, and a removable grid.

This x-ray solution can fit in virtually any x-ray room and drastically cuts down installation time reducing clinical downtime.



The Acuity CDR software is fully integrated to the x-ray generator and x-ray techniques are seamlessly preselected based on the anatomical region of the exam.

Acuity CDR



## BUCKY ENCLOSURE



Features fully counterbalanced powerful magnetic lock brakes from the handle control to raise and lower the system effortlessly. The bucky also rotates to maintain tube alignment during arm angulation positioning.



## TUBE HEAD



Features a dual laser rotating LED collimator with precise alignment and a fixed 40in. SID system with optional Automatic Exposure Control (AEC).



## HEAD TO TOE POSITIONING



The Acuity CDR can seamlessly raise and lower to accommodate all standing exams from skull to lower extremities. The arm utilizes 20° angulation which offsets anatomical anomalies for higher detailed x-rays and better patient care.

