



CFPM300 & CFPM301 Imaging Tables

Oakworks Imaging Tables have a large height range helping to create safe patient transfers while allowing doctors to stand comfortably. Our large lifting column is very rigid, ensuring rock solid support for fast imaging and procedures.



Integrated Headrest Table

Motions

- Motion 1 - Height Range 26" - 43" (CFPM300)
- Motion 1 - Height Range 25" - 43" (CFPM301)
- Motion 2 - Trendelenburg Tilt +15/-11°
- Motion 3 - Longitudinal Travel 10" (CFPM301)
- Motion 3 - Lateral Tilt $\pm 15^\circ$ (CFPM300)

Features

- Choice of integrated headrest or rectangular carbon fiber tops with 56" of metal free imaging
- 2" Radiolucent Comfort Foam[®] table top pad
- Warranty is 4 years parts and 3 years labor

SPECIFICATIONS

- Large Adjustable Height Range 26" - 43" (CFPM300), 25" - 43" (CFPM301)
- 56" Metal Free Imaging Space
- Cost Effective 3 Movement Table
- 500 lb Patient Weight Capacity
- Warranty: 4 years parts, 3 years labor

CFPM300 #81803 Integrated #81810 Rectangular CFPM301 #81804 Integrated #81811 Rectangular

INCLUDED FEATURES	
Carbon Fiber Top 56" Metal free imaging space	Integrated Headrest 22" (56cm) wide x 84" (214cm) long Rectangular 24" (61cm) wide x 80" (203cm) long
CFPM300 MOVEMENTS	
Motion 1	Electric Height Range 26" - 43" (66cm-109cm) (without pad)
Motion 2	Electric Trendelenburg tilt +15°/-11°
Motion 3	Electric Lateral tilt ±15°
CFPM301 MOVEMENTS	
Motion 1	Electric Height Range 25" - 43" (64-109cm) (without pad)
Motion 2	Electric Trendelenburg tilt +15°/-11°
Motion 3	Electric Longitudinal travel 10" (25cm)
Controls	All functions hand control
Weight Capacity	500 lb (227kg)
Product Weight	475 lb (216kg)
Included Features	2" (5cm) Radiolucent Comfort Foam™ table top pad, 4" (10cm) locking casters, Crescent face rest pad (included with integrated headrest top), patient restraint straps (2)
Warranty	4 years parts, 3 years labor
Safety Listings	FDA listed, CE marked, cETLus listed, CB Certificate

OPTIONS *must be ordered with table.*

5" (13cm) casters	Upgrade from the standard 4" (10cm)
Foot Control	All functions control
Non-USA electrical systems	Please inquire

ACCESSORIES *order any time.*

		#
Arm Hammock	Use with Integrated Headrest Top	48751
Carbon Fiber Arm Board (6" wide)	180° Articulation	83575
Fluoro Extender	Polycarbonate	33975
Interventional Supine Pillow System	Contoured pad and support wedge pad	81821
Spine Positioning System II	8 piece system	81819
T-rail Clamp	6" (15cm) long rail	63336



MEDICAL X-RAY SYSTEMS
EQUIPMENT • SUPPLIES • SERVICE
800-356-3388

Associated X-Ray Imaging Corp.
49 Newark Street
Haverhill, MA 01832
Toll Free: 800-356-3388
Phone: 978-374-6371
Fax: 978-521-2214
www.associatedxray.com
sales@associatedxray.com



Rectangular Top Table

ACCESSORIES

SPINE POSITIONING SYSTEM II

Radiolucency & adjustability for ideal cervical, thoracic & lumbar visualization. Minimize unwanted movement & reduce procedure time.



ARM HAMMOCK

This unique product will comfortably & easily support your patient's arms when prone. Easily attached to the Integrated Head Rest table top.



SINGLE CARBON FIBER ARM BOARD

Fully radiolucent with a 180° range of motion, the arm board requires no mounting or attachment hardware and simply slides under the table top pad. It can be used on any fluoroscopy table & will support arms and feet in a variety of positioning configurations.



FLUORO EXTENDER

Made of radiolucent, polycarbonate material, it's organic shape is conducive to a variety of table placement applications. 30" long x 23.5" wide x 1/4" thick



6" ADJUSTABLE T-RAIL CLAMP

The adjustable T-rail clamp can be positioned along the grooved table top edge. Clamp is also removable.



5" CASTER (OPTIONAL)

Recommended for mobile service providers.

