



SPECIFICATIONS

Weight		
With fixed holder	35kg (77lbs)	
With tiltable holder	38kg (84lbs)	
Max load on detector holder		
	25kg (55lbs)	
	<i>Fixed</i>	<i>Tilt</i>
Max height (top of active image surface)	1810mm (71,3")	1900mm (74,8")
Min height (bottom of active image surface)	56mm (2,2")	18,5mm (0,7")
Stroke:	1324mm (52,1")	1451,5mm (57,1")
Detector sizes	35 x 43cm (14" x 17"), 43 x 43cm (17" x 17")	
Max speed:	125mm/sec	
Input power:	100-240V AC, 50-60Hz	
Battery:	24V, 7.2h, Charging time: 4 hours (80%) 5 hours (100%) Operation time: More than a week	
IP class:	IP41	
Certifications:	CE, MD 2006/42/EG IIA	

OPTIONS

Fixed or tiltable detector holder
Lateral patient support handle
Auto-tracking w/tube head(require enabled tracking function on M1)
Docking for transport, available when used w/M1

! Solutions for tomorrow

Find out more about **TRUE X-RAY MOBILITY®** at www.solutionsfortomorrow.com www.mobilewallstand.com

Associated X-Ray Imaging Corp.
 49 Newark Street,
 Haverhill MA 01832
 800-356-3388
 sales@associatedxray.com
 associatedxray.com

MEDICAL X-RAY SYSTEMS
 EQUIPMENT • SUPPLIES • SERVICE

Solutions for tomorrow AB reserves the right to change specifications without prior notice.

© 2022 All right reserved.

Solutions for tomorrow AB, Väckelsång, Sweden

Developed and manufactured in Sweden

! Solutions for tomorrow

SETTING NEW STANDARDS WITH

MOBILE WALL STAND



Product_flyer_!Mobile_wall_stand_2p_A1_EN



!MOBILE WALL STAND

Solutions for tomorrow proudly present our unique motorized and battery operated !Mobile wall stand.

MOBILE AND MOTORIZED

- Small and light
- Easy to move
- Battery operated
- Height of the detector is adjusted with wired hand control
- Wireless auto-tracking of the tube head when used with SFT M1

FLEXIBLE

- Can be used with any equipment
- Fits all types of detectors and cassettes
- Use with or without grid
- Long stroke (up/down reach)
- Covers all RAD examinations (with tilt option)

COST-EFFICIENT

- Turnkey solution
- Save time and cost of installation, no installation needed
- Ready to use within seconds



READY TO BE USED WITH ANY
PORTABLE X-RAY BRAND