

# OPAL-VET

Same great enterprise-level PACS features with customized functionalities for the veterinary practice. Add Viztek's DR for complete digital workflow.

- Viztek's PACS is Web-based and can be accessed securely from anywhere over an internet connection
- Import CDs from outside facilities
- Exam room viewing software allows the freedom to view in multiple exam rooms
- Multi-modality, multi-vendor compatibility
- CD Burning, DICOM routing and DICOM print functions
- Scan documents into exams
- Simple, fast, robust
- Pet name and owner tracking
- Fields for species and breeds
- Field for spayed or neutered status



2217 US Highway 70E • Garner, NC 27529  
Phone: 800.366.5343 • Fax: 919.772.2530  
Email: [info@viztek.net](mailto:info@viztek.net) • [www.viztek.net](http://www.viztek.net)

Viztek and Opal-VET are trademarks of Viztek LLC. © 2010 Viztek LLC All rights reserved.

## DR Features

### Flat Panel Detector Assembly

- Express DR
- Detector Housing
- 17"x17" imaging area
- 3.6 lp/mm
- Off-center imaging
- 9 million pixels



### Panel Features

- Unsurpassed Flat Panel Technology
- Superior Image Quality
- High Sensitivity
- Sleek Slim Design

### Panel Specs

- Pixel Size - 143
- Pixel Array - 3072 x 3072
- Pixel Area - 17in x 17in

### CR Features

- Single Cassette Feed
- 4 Cassettes Included
- Opal Front-End Viewing Software
- Data Acquisition : 20 bits/pixel
- Pixel Matrix Size : 3480x4248
- 14 x 17 cassette size : 10 pixels/mm
- Power Requirement : 120V/60Hz
- Small Table-Top Footprint
- Dimension : 27.2 x 27.6 x 18.2"



**MEDICAL X-RAY SYSTEMS**  
EQUIPMENT • SUPPLIES • SERVICE

800-356-3388  
978-374-6371  
Fax - 978-521-2214

49 Newark Street  
Haverhill, MA 01832  
[sales@associatedxray.com](mailto:sales@associatedxray.com)  
[www.associatedxray.com](http://www.associatedxray.com)