



by **SEDECAL**



Premium Vet *The premium choice in veterinary specific x-ray.*

High Frequency Compact Generator • Anatomical Programming (APR) • Variable Focal Distance
Moveable Grid Cabinet • 4 Way Float Top Table • Heavy Gauge Steel Construction
Rotor Brake • Capable of off table exposures • Horizontal Beam Imaging

High Frequency X-ray Generator

400 mA
125 kVp
500mAs
Serial interface control
Push button console with APR window
Anatomical Programmed Radiography (APR)
208/220/240 VAC Single Phase Power
Exposure Foot pedal

Variable Tube Stand and Table

Variable SID
Tube stand independent of table
Rotating tube head and collimator
5 foot long 4 way float top table
Foot treadle to operate table locks
Hands free collimator light control
Manual Collimator with light switch
Grid Cabinet with 8:1 103LP Grid
Stainless steel cassette tray
Longitudinal movement of tube stand
Wall cassette holder for off table exams

X-ray Tube – Rotating Anode

1.0/2.0 Focal Spots
16° Target Angle
140,000 HU
6 meter HV cables

Optional Features

6 foot Table Top
Fractional Focal Spot X-ray tube
Non-certified medical collimator
3 phase and battery powered generators

800.920.9525
Info@VetRay.com

www.VetRay.com



