



Your Practice. Your Preference.
Your Productivity Solution.

LEGGERA DR SOLUTIONS

Leggera DR Solutions are available in a cesium wireless cassette size wifi panel and a 4143 fixed panel. Both offer superior image quality, low dose, and offer images within seconds of exposure for speedy throughput.

The 35x43cm wireless panel features an ultra-light design and is highly versatile, fitting in various x-ray room environments. At only 6.17lbs and with up-to 8-hours of battery life, the (ISO 4090) cassette size panel uses ethernet communication technology and includes a back-up connection cable. An optional battery charger is available and can charge up to three batteries at a time.

The fixed 4143 panel is an ideal solution for various x-ray room configurations including the FIT Kit. This incorporates the Fixed and Wireless panel in a dual room.

Leggera DR Solutions	FIXED	WIRELESS
Pixel Size	148um	148um
Dimensions	19.7x19.3x1.8 in	ISO 4090 Cassette Size
Weight	27.55lbs	6.17lbs
DQE @ 0lp/mm	64%	66%
AD Conversion	16 Bits	16 Bits
Image Display Time	Less than 5 sec	Preview in 3 sec
X-Ray Gen Voltage Range	40-150kv	40-150kv



2217 US HWY 70 East • Garner, NC 27529 • Phone: 800.366.5343
Email: info@viztek.net • www.viztek.net

Viztek and Leggera DR are trademarks of Viztek LLC. © 2012 Viztek LLC. All rights reserved.



800-356-3388
978-374-6371
Fax – 978-521-2214

49 Newark Street
Haverhill, MA 01832
sales@associatedxray.com
www.associatedxray.com

