



Xmaru VetView Veterinary Console

VETERINARY HOSPITALS • VETERINARY MOBILE UNITS • VETERINARY CLINICS • VETERINARY SCHOOLS



XmaruVetView Console

Rayence introduces their latest acquisition software for veterinary, the Xmaru VetView. The Xmaru VetView console is at the heart of what drives the Xmaru Panel Series to its full potential. Xmaru VetView uses the industry's best animal image algorithm parameters for the highest quality diagnostic images. Designed with workflow and efficiency in mind, the Xmaru VetView was developed specifically for the veterinary environments.

The Xmaru VetView comes with standard features like DICOM Modality Work List (MWL), DICOM Printing, and DVD/CD Burning with viewer. We provide you the ability export images as JPEG, TIF, GIF, & BITMAP. The Xmaru VetView comes with image mirroring, backup and restoration solutions for disaster recovery on our workstations.

Xmaru VetView is today's leading veterinary console software. Xmaru VetView comes with a five year warranty with free updates and upgrades. With a simple work flow and high end image quality, the XmaruVetView software is the best in its class for today's veterinarians

Xmaru VetView Console Standard Features

DICOM MWL	JPEG Import	Left & Right Marker	DICOM Send
Exposure Index	DICOM Import	Shuttering	Query & Retrieve
Image Rotation	Zoom	Auto W/L ROI	CD/DVD Burning with Viewer
Flip Vertical	Panning	Image Reprocessing	Backup & Restore Function
Flip Horizontal	Annotations	DICOM Export	Mirror Backup
1 to 1 True Size	Predefine Annotations	JPEG Export	Hanging Protocol
Magnification	Line Measurements	Arrow Pointers	DICOM Printing
Veterabab Heart Score	Clock Face VD	Clock Face Lateral	Notes

XmaruVetView Console

Available Languages

- English
- Spanish
- Japanese
- Chinese - Simplified
- Chinese - Traditional
- French
- German
- Russian
- Thai



800-356-3388

978-374-6371

Fax - 978-521-2214

49 Newark Street
 Haverhill, MA 01832

sales@associatedxray.com

www.associatedxray.com



Rayence Inc.

440 Sylvan Avenue, Suite 220

Englewood Cliffs, NJ 07632

P: (201) 569-0445

F: (201) 569-0448

www.rayenceusa.com